

Initial F Lapidus™: single use sterile kits for plantar lapidus arthrodesis

Initial F Lapidus™ - Kits references




Indications:

The implants of the Initial F Lapidus range™ are intended for arthrodeses, fractures and osteotomies fixation and revision surgeries of the foot in adults.

Contraindications:

- Serious vascular deterioration, bone devitalization,
- Pregnancy.
- Acute or chronic local or systemic infections.
- Lack of musculo-cutaneous cover, severe vascular deficiency affecting the concerned area.
- Insufficient bone quality preventing a good fixation of the implants into the bone,
- Muscular deficit, neurological deficiency or behavioral disorders, which could submit the implant to abnormal mechanical strains.
- Allergy to one of the materials used or sensitivity to foreign bodies.
- Serious problems of non-compliance, mental or neurological disorders, failure to follow post-operative care recommendations.
- Unstable physical and/or mental condition.



NEWCLIP-TECHNICS



Initial F Lapidus™: single use sterile kits for plantar lapidus arthrodesis

Initial F Lapidus™ - Kits references

> Kits content - Instruments

STERILE INSTRUMENTATION FOR PLANTAR LAPIDUS PLATE

Description	Content	Qty
Lapidus Arthrodesis kit - Plantar plate - Size 1	5.8 mm single use handle	1
	Handle for guide gauge	1
	Length gauge	1
	Ø2.7 mm threaded guide gauge for Ø3.5 mm screws	1
	Ø2.7 mm quick coupling drill bit - L 125 mm	1
	T8 quick coupling screwdriver	1
	Pin Ø1.6 L100 mm	1

STERILE INSTRUMENTATION FOR Ø4.0 MM CANNULATED SCREWS*

Description	Content	Qty
Kit for Ø4.0 mm cannulated headed screws	Length gauge for pin Ø1.3 mm	1
	2 in 1: 2.5 mm hexagonal screwdriver - Ø6.0 mm countersink	1
	Ø2.9 mm drill bit - cannula 1.4 mm - L 120 mm - AO Ø4.5 mm quick coupling	1
	Pin Ø1.3 L120 mm	3



* The Ø4.0 mm cannulated headed screws (Ref. H1.4QT4.0xxL-STI) are available in supplemental sterile screw caddy

> Kits content - Implants

KIT-LV1D-CQ4.0 or KIT-LV1G-CQ4.0 Lapidus Arthrodesis kit - Left or Right - Size 1

Ref.	Description	KIT-LV1D-CQ4.0 or KIT-LV1G-CQ4.0
FLTxDV1	Plantar Lapidus Arthrodesis plate - Narrow - Size 1 - Right or Left	1
RLT3.5L20	Non locking screw - Ø3.5mm - L 20 mm	1
RLT3.5L26	Non locking screw - Ø3.5 mm - L 26 mm	1
RLT3.5L28	Non locking screw - Ø3.5 mm - L 28 mm	1
SLT3.5L14	Locking screw - Ø3.5 mm - L14 mm	1
SLT3.5L16	Locking screw - Ø3.5 mm - L16 mm	1
SLT3.5L20	Locking screw - Ø3.5 mm - L20mm	1
SLT3.5L26	Locking screw - Ø3.5 mm - L26mm	1

* Each kit includes one plate (left or right).

NB : Additional sterile screws can be ordered separated (Supplemental sterile screw caddy).

> Additional instrumentation kits

INITIAL F™ - LAPIDUS ARTHRODESIS TEMPLATES

Ref.	Description
ANC958	FLTDVx and FLTGVx templates

> Removal and Rescue kits

INITIAL F™ - REMOVAL AND RESCUE KITS

Ref.	Description
KIT-RESCUE-6	Rescue kit for Initial F lapidus
KIT-REMOVE-2	Removal kit for T8 hexalobe



NEWCLIP TECHNICS
PA de la Lande Saint Martin - 45 rue des Garottières
44115 Haute Goulaine (France)
Tél. : +33 (0)2 28 21 37 12 - Fax : +33 (0)2 40 63 68 37
Mail : orders@newcliptechnics.com
www.newcliptechnics.com

NEWCLIP USA
642 Larkfield Center
Santa Rosa CA 95403, USA
Phone : + 1 707 230 5078
Mail : customerservice@newclipusa.com
www.newclipusa.com

The information presented in this brochure is intended to demonstrate a Newclip Technics product. Always refer to the package insert, product label and/or user instructions before using any Newclip Technics product. Surgeons must always rely on their own clinical judgment when deciding which products and techniques to use with their patients. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Newclip Technics representative if you have questions about the availability of Newclip Technics products in your area.

